

HUMAN LAMP3 PROTEIN, HIS TAG

目录: 11807

产品名称: Human LAMP3 Protein

规格: 10 µg, 50 µg and 100 µg

基因符号: CD208;DC-LAMP;DC LAMP;DCLAMP;LAMP;LAMP-3;TSC403

Target: LAMP3

UNIPROT ID: Q9UQV4

描述: Recombinant human LAMP3 protein with C-terminal 6xHis tag

背景: Dendritic cells (DCs) are the most potent antigen-presenting cells. Immature DCs efficiently capture antigens and differentiate into interdigitating dendritic cells (IDCs) in lymphoid tissues that induce primary T-cell responses (summary by de Saint-Vis et al., 1998 [PubMed 9768752]).[supplied by OMIM, Dec 2010]

Species/Host: HEK293

Molecular Weight: The protein has a predicted molecular mass of 38.6 kDa after removal of the signal peptide. The apparent molecular mass of LAMP3-His is approximately 35-130 kDa due to glycosylation.

Molecular Characterization: LAMP3(Lys28-Thr381) 6×His tag

纯化: The purity of the protein is greater than 95% as determined by SDS-PAGE and Coomassie blue staining.

Formulation & Reconstitution: Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% – 8% trehalose is added as protectants before lyophilization.

储存和运输: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.



Figure 1. Human LAMP3 Protein, His Tag on SDS-PAGE under reducing condition.