

HUMAN LDLR PROTEIN, HIS TAG

目录: 11859

产品名称: Human LDLR Protein

规格: 10 µg, 50 µg and 100 µg

基因符号: FH;FHC;FHCL1;LDLCQ2

Target: LDLR

UNIPROT ID: P01130

描述: Recombinant Human LDLR Protein with C-terminal 6xHis tag

背景: The low density lipoprotein receptor (LDLR) gene family consists of cell surface proteins involved in receptor-mediated endocytosis of specific ligands. Low density lipoprotein (LDL) is normally bound at the cell membrane and taken into the cell ending up in lysosomes where the protein is degraded and the cholesterol is made available for repression of microsomal enzyme 3-hydroxy-3-methylglutaryl coenzyme A (HMG CoA) reductase, the rate-limiting step in cholesterol synthesis. At the same time, a reciprocal stimulation of cholesterol ester synthesis takes place. Mutations in this gene cause the autosomal dominant disorder, familial hypercholesterolemia. Alternate splicing results in multiple transcript variants.[provided by RefSeq, Sep 2010]

Species/Host: HEK293

Molecular Weight: The protein has a predicted molecular mass of 85.6 kDa after removal of the signal peptide. The apparent molecular mass of LDLR-His is approximately 100-130 kDa due to glycosylation.

Molecular Characterization: LDLR(Ala22-Arg788) 6xHis tag

纯化:: The purity of the protein is greater than 85% as determined by SDS-PAGE and Coomassie blue staining.

Formulation & Reconstitution: Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% – 8% trehalose is added as protectants before lyophilization.

储存和运输: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.




Figure 1. Human LDLR Protein, His Tag on SDS-PAGE under reducing condition.