

HUMAN LYPD3(31-286) PROTEIN, HIS TAG

Cat.#: 11890

Product Name: Human LYPD3(31-286) Protein

Size: 10 µg, 50 µg and 100 µg

Synonyms: MIG-C4;4.4A

Target: LYPD3

UNIPROT ID: O95274

Description: Recombinant Human LYPD3(31-286) Protein with C-terminal 6xHis tag

Background: Supports cell migration. May be involved in urothelial cell-matrix interactions. May be involved in tumor progression.[UniProtKB/Swiss-Prot Function]

Species/Host: HEK293

Molecular Weight: The protein has a predicted molecular mass of 27.7 kDa after removal of the signal peptide. The apparent molecular mass of LYPD3(31-286)-His is approximately 35-70 kDa due to glycosylation.

Molecular Characterization: LYPD3(Leu31-His286) 6×His tag

Purity: The purity of the protein is greater than 85% as determined by SDS-PAGE and Coomassie blue staining.

Formulation & Reconstitution: Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% – 8% trehalose is added as protectants before lyophilization.

Storage & Shipping: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.



Figure 1. Human LYPD3(31-286) Protein, His Tag on SDS-PAGE under reducing condition.