

## HUMAN M-CSF PROTEIN, HIS TAG

目录: 11369

产品名称: Human M-CSF Protein

规格: 10 µg, 50 µg and 100 µg

基因符号: CSF-1;MCSF

**Target:** M-CSF

**UNIPROT ID:** P09603

**描述:** Recombinant Human M-CSF with C-terminal 6xHis tag

**背景:** The protein encoded by this gene is a cytokine that controls the production, differentiation, and function of macrophages. The active form of the protein is found extracellularly as a disulfide-linked homodimer, and is thought to be produced by proteolytic cleavage of membrane-bound precursors. The encoded protein may be involved in development of the placenta. Alternate splicing results in multiple transcript variants. [provided by RefSeq, Sep 2011]

**Species/Host:** HEK293

**Molecular Weight:** The protein has a predicted molecular mass of 26 kDa after removal of the signal peptide. The apparent molecular mass of M-CSF-His is approximately 25-55 kDa due to glycosylation.

**Molecular Characterization:** M-CSF(Glu33-Arg255) 6×His tag

**纯化:** The purity of the protein is greater than 95% as determined by SDS-PAGE and Coomassie blue staining.

**Formulation & Reconstitution:** Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% – 8% trehalose is added as protectants before lyophilization.

**储存和运输:** Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.




Figure 1. Human M-CSF Protein, His Tag on SDS-PAGE under reducing condition.

## Human M-CSF, His Tagged protein ELISA

0.2  $\mu$ g of Human M-CSF, His tagged protein per well




Figure 2. ELISA plate pre-coated by 2  $\mu$ g/ml (100  $\mu$ l/well) Human M-CSF Protein, His Tag (11369) can bind M-CSF Neutralizing antibody 28104 in a linear range of 2.44-156.25 ng/mL.




Figure 3. Flow cytometry analysis with 1  $\mu$ g/ml Human M-CSF Protein, His tag (11369) on Expi293 cells transfected with human CSFIR (Blue histogram) or Expi293 transfected with irrelevant protein (Red histogram).