

## HUMAN MERTK PROTEIN, HIS TAG

**Cat.#:** 11262

**Product Name:** Human MERTK Protein

**Size:** 10 µg, 50 µg and 100 µg

**Synonyms:** MER;RP38;c-Eyk;c-mer;Tyro12

**Target:** MERTK

**UNIPROT ID:** Q12866

**Description:** Recombinant human MERTK Protein with C-terminal 10xHis tag

**Background:** This gene is a member of the MER/AXL/TYRO3 receptor kinase family and encodes a transmembrane protein with two fibronectin type-III domains, two Ig-like C2-type (immunoglobulin-like) domains, and one tyrosine kinase domain. Mutations in this gene have been associated with disruption of the retinal pigment epithelium (RPE) phagocytosis pathway and onset of autosomal recessive retinitis pigmentosa (RP). [provided by RefSeq, Jul 2008]

**Species/Host:** HEK293

**Molecular Weight:** The protein has a predicted molecular mass of 54.0 kDa after removal of the signal peptide.

**Molecular Characterization:** MERTK(Ala21-Ile505) 10×His tag

**Purity:** The purity of the protein is greater than 85% as determined by SDS-PAGE and Coomassie blue staining.

**Formulation & Reconstitution:** Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% – 8% trehalose is added as protectants before lyophilization.

**Storage & Shipping:** Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.




Figure 1. Human MERTK Protein, His Tag on SDS-PAGE under reducing condition.