

**HUMAN MICAA3 PROTEIN****Cat.#:** 12283**Product Name:** Human MICAA3 Protein**Size:** 10 µg, 50 µg and 100 µg**Synonyms:** MICA;MIC-A;PERB11.1**Target:** MICAA3**UNIPROT ID:** Q29983**Description:** Recombinant human MICAA3 Protein with C-terminal 6xHis tag**Background:** This gene encodes the highly polymorphic major histocompatibility complex class I chain-related protein A. The protein product is expressed on the cell surface, although unlike canonical class I molecules it does not seem to associate with beta-2-microglobulin. It is a ligand for the NKG2-D type II integral membrane protein receptor. The protein functions as a stress-induced antigen that is broadly recognized by intestinal epithelial gamma delta T cells. Variations in this gene have been associated with susceptibility to psoriasis I and psoriatic arthritis, and the shedding of MICA-related antibodies and ligands is involved in the progression from monoclonal gammopathy of undetermined significance to multiple myeloma. Alternative splicing of this gene results in multiple transcript variants. [provided by RefSeq, Jan 2014]**Species/Host:** HEK293**Molecular Weight:** The protein has a predicted molecular mass of 12.5 kDa after removal of the signal peptide.**Molecular Characterization:** MICAA3(Arg203-His306) 6xHis tag**Purity:** The purity of the protein is greater than 85% as determined by SDS-PAGE and Coomassie blue staining.**Formulation & Reconstitution:** Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% – 8% trehalose is added as protectants before lyophilization.**Storage & Shipping:** Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.




Figure 1. Human MICA±3 Protein, His Tag on SDS-PAGE under reducing condition.