

HUMAN MMP11 PROTEIN, HFC TAG

Cat.#: 11703

Product Name: Human MMP11 Protein

Size: 10 µg, 50 µg and 100 µg

Synonyms: MMP-11;SL-3;ST3;STMY3

Target: MMP11

UNIPROT ID: P24347

Description: Recombinant human MMP11 protein with C-terminal human Fc tag

Background: Proteins of the matrix metalloproteinase (MMP) family are involved in the breakdown of extracellular matrix in normal physiological processes, such as embryonic development, reproduction, and tissue remodeling, as well as in disease processes, such as arthritis and metastasis. Most MMP's are secreted as inactive proproteins which are activated when cleaved by extracellular proteinases. However, the enzyme encoded by this gene is activated intracellularly by furin within the constitutive secretory pathway. Also in contrast to other MMP's, this enzyme cleaves alpha 1-proteinase inhibitor but weakly degrades structural proteins of the extracellular matrix. [provided by RefSeq, Jul 2008]

Species/Host: HEK293

Molecular Weight: The protein has a predicted molecular mass of 76.5 kDa after removal of the signal peptide. The apparent molecular mass of MMP11-hFc is approximately 70-100 kDa due to glycosylation.

Molecular Characterization: MMP11(Arg32-Leu488) hFc(Glu99-Ala330)

Purity: The purity of the protein is greater than 95% as determined by SDS-PAGE and Coomassie blue staining.

Formulation & Reconstitution: Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% – 8% trehalose is added as protectants before lyophilization.

Storage & Shipping: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.




Figure 1. Human MMP11 Protein, hFc Tag on SDS-PAGE under reducing condition.