

HUMAN MSMP PROTEIN, HFC TAG

目录: 11837

产品名称: Human MSMP Protein

规格: 10 µg, 50 µg and 100 µg

基因符号: PSMP

Target: MSMP

UNIPROT ID: Q1L6U9

描述: Recombinant Human MSMP Protein with C-terminal human Fc tag

背景: This gene encodes a member of the beta-microseminoprotein family. Members of this protein family contain ten conserved cysteine residues that form intra-molecular disulfide bonds. The encoded protein may play a role in prostate cancer tumorigenesis. [provided by RefSeq, Jan 2011]

Species/Host: HEK293

Molecular Weight: The protein has a predicted molecular mass of 37.3 kDa after removal of the signal peptide. The apparent molecular mass of MSMP-hFc is approximately 35–55 kDa due to glycosylation.

Molecular Characterization: MSMP(Ser36–Ser139) hFc(Glu99–Ala330)

纯化:: The purity of the protein is greater than 95% as determined by SDS-PAGE and Coomassie blue staining.

Formulation & Reconstitution: Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% – 8% trehalose is added as protectants before lyophilization.

储存和运输: Store at –20°C to –80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at –80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.



Figure 1. Human MSMP Protein, hFc Tag on SDS-PAGE under reducing condition.