

HUMAN MXRA8 PROTEIN, HIS TAG

目录: 11984

产品名称: Human MXRA8 Protein

规格: 10 µg, 50 µg and 100 µg

基因符号: ASP3

Target: MXRA8

UNIPROT ID: Q9BRK3

描述: Recombinant Human MXRA8 Protein with C-terminal 6xHis tag

背景: Predicted to be involved in establishment of glial blood-brain barrier. Located in extracellular exosome. [provided by Alliance of Genome Resources, Apr 2022]

Species/Host: HEK293

Molecular Weight: The protein has a predicted molecular mass of 36.3 kDa after removal of the signal peptide. The apparent molecular mass of MXRA8-His is approximately 35–55 kDa due to glycosylation.

Molecular Characterization: MXRA8(Val20–Gln340) 6×His tag

纯化:: The purity of the protein is greater than 85% as determined by SDS-PAGE and Coomassie blue staining.

Formulation & Reconstitution: Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% – 8% trehalose is added as protectants before lyophilization.

储存和运输: Store at –20°C to –80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at –80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.




Figure 1. Human MXRA8 Protein, His Tag on SDS-PAGE under reducing condition.