

HUMAN MYELOBLASTIN PROTEIN, HFC TAG

目录: 11706

产品名称: Human Myeloblastin Protein

规格: 10 µg, 50 µg and 100 µg

基因符号: ACPA;AGP7;C-ANCA;CANCA;MBN;MBT;NP-4;NP4;P29;PR-3;PR3

Target: Myeloblastin

UNIPROT ID: P24158

描述: Recombinant human Myeloblastin protein with C-terminal human Fc tag

背景: Polymorphonuclear leukocyte serine protease that degrades elastin, fibronectin, laminin, vitronectin, and collagen types I, III, and IV (in vitro) and causes emphysema when administered by tracheal insufflation to hamsters. [UniProtKB/Swiss-Prot Function]

Species/Host: HEK293

Molecular Weight: The protein has a predicted molecular mass of 49.8 kDa after removal of the signal peptide.

Molecular Characterization: Myeloblastin(Ile28-Arg248) hFc(Glu99-Ala330)

纯化: The purity of the protein is greater than 95% as determined by SDS-PAGE and Coomassie blue staining.

Formulation & Reconstitution: Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% – 8% trehalose is added as protectants before lyophilization.

储存和运输: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.