

HUMAN NINJ1 PROTEIN, HFC TAG

目录: 11906

产品名称: Human NINJ1 Protein

规格: 10 µg, 50 µg and 100 µg

基因符号: NIN1;NINJURIN

Target: NINJ1

UNIPROT ID: Q92982

描述: Recombinant Human NINJ1 Protein with C-terminal human Fc tag

背景: The ninjurin protein is upregulated after nerve injury both in dorsal root ganglion neurons and in Schwann cells (Araki and Milbrandt, 1996 [PubMed 8780658]). It demonstrates properties of a homophilic adhesion molecule and promotes neurite outgrowth from primary cultured dorsal root ganglion neurons.[supplied by OMIM, Aug 2009]

Species/Host: HEK293

Molecular Weight: The protein has a predicted molecular mass of 34.6 kDa after removal of the signal peptide. The apparent molecular mass of NINJ1-hFc is approximately 35-55 kDa due to glycosylation.

Molecular Characterization: NINJ1(Met1-Leu80) hFc(Glu99-Ala330)

纯化: The purity of the protein is greater than 95% as determined by SDS-PAGE and Coomassie blue staining.

Formulation & Reconstitution: Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% – 8% trehalose is added as protectants before lyophilization.

储存和运输: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.



Figure 1. Human NINJ1 Protein, hFc Tag on SDS-PAGE under reducing condition.