

HUMAN OSMR PROTEIN, HIS TAG

目录: 11782

产品名称: Human OSMR Protein

规格: 10 µg, 50 µg and 100 µg

基因符号: IL-31R-beta;IL-31RB;OSMRB;OSMRbeta;PLCA1

Target: OSMR

UNIPROT ID: Q99650

描述: Recombinant human OSMR protein with C-terminal 6xHis tag

背景: This gene encodes a member of the type I cytokine receptor family. The encoded protein heterodimerizes with interleukin 6 signal transducer to form the type II oncostatin M receptor and with interleukin 31 receptor A to form the interleukin 31 receptor, and thus transduces oncostatin M and interleukin 31 induced signaling events. Mutations in this gene have been associated with familial primary localized cutaneous amyloidosis. Alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Dec 2009]

Species/Host: HEK293

Molecular Weight: The protein has a predicted molecular mass of 81.9 kDa after removal of the signal peptide. The apparent molecular mass of OSMR-His is approximately 100-250 kDa due to glycosylation.

Molecular Characterization: OSMR(Glu28-Met740) 6×His tag

纯化:: The purity of the protein is greater than 85% as determined by SDS-PAGE and Coomassie blue staining.

Formulation & Reconstitution: Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% – 8% trehalose is added as protectants before lyophilization.

储存和运输: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.




Figure 1. Human OSMR Protein, His Tag on SDS-PAGE under reducing condition.