

HUMAN PDGF-BB PROTEIN

Cat.#: 12035

Product Name: Human PDGF-BB Protein

Size: 10 µg, 50 µg and 100 µg

Synonyms: PDGFBB;PDGF-BB

Target: PDGF-BB

UNIPROT ID: P01127

Description: Recombinant Human Platelet-Derived Growth Factor BB is produced by our E.coli expression system and the target gene encoding Ser82-Thr190 is expressed.

Background: Platelet-Derived Growth Factor Subunit B (PDGFB) belongs to the PDGF/VEGF growth factor family. Platelet-derived growth factor is a potent mitogen for cells of mesenchymal origin. PDGFB can exist either as a homodimer (PDGF-BB) or as a heterodimer with the platelet-derived growth factor alpha polypeptide (PDGF-AB), where the dimers are connected by disulfide bonds. Mutations in this gene are associated with meningioma. Binding of PDGFB to its receptor elicits a variety of cellular responses. In addition, PDGFB is released by platelets upon wounding and plays an important role in stimulating adjacent cells to grow and thereby heals the wound.

Species/Host: E.coli

Molecular Weight: 12.42 KDa

Molecular Characterization: Not available

Purity: Greater than 98% as determined by reducing SDS-PAGE.

Formulation & Reconstitution: Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% – 8% trehalose is added as protectants before lyophilization.

Storage & Shipping: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.




Figure 1. Greater than 95% as determined by reducing SDS-PAGE.