

HUMAN PLA2R1(21-164+223-359+1097-1246) PROTEIN

Cat.#: 12256

Product Name: Human PLA2R1(21-164+223-359+1097-1246) Protein

Size: 10 µg, 50 µg and 100 µg

Synonyms: CLEC13C;PLA2-R;PLA2G1R;PLA2IR;PLA2R

Target: PLA2R1

UNIPROT ID: Q13018

Description: Recombinant Human PLA2R1(21-164 223-359 1097-1246) Protein with C-terminal 10xHis tag

Background: This gene represents a phospholipase A2 receptor. The encoded protein likely exists as both a transmembrane form and a soluble form. The transmembrane receptor may play a role in clearance of phospholipase A2, thereby inhibiting its action. Polymorphisms at this locus have been associated with susceptibility to idiopathic membranous nephropathy. Alternatively spliced transcript variants encoding different isoforms have been identified.[provided by RefSeq, Sep 2010]

Species/Host: HEK293

Molecular Weight: The protein has a predicted molecular mass of 50.3 kDa after removal of the signal peptide. The apparent molecular mass of PLA2R1(21-164 223-359 1097-1246)-10xHis is approximately 35-70 kDa due to glycosylation.

Molecular Characterization: PLA2R1(Glu21-Lys164 Thr223-Asn359 Glu1097-Leu1246) 10xHis tag

Purity: The purity of the protein is greater than 85% as determined by SDS-PAGE and Coomassie blue staining.

Formulation & Reconstitution: Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% – 8% trehalose is added as protectants before lyophilization.

Storage & Shipping: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.



Figure 1. Human PLA2R1(21-164 223-359 1097-1246) Protein, His Tag on SDS-PAGE under reducing condition.