

## HUMAN PRTN3 PROTEIN, HIS TAG

目录: 11817

产品名称: Human PRTN3 Protein

规格: 10 µg, 50 µg and 100 µg

基因符号: ACPA;AGP7;C-ANCA;CANCA;MBN;MBT;NP-4;NP4;P29;PR-3;PR3

**Target:** PRTN3

**UNIPROT ID:** P24158

**描述:** Recombinant Human PRTN3 Protein with C-terminal 6xHis tag

**背景:** Serine protease that degrades elastin, fibronectin, laminin, vitronectin, and collagen types I, III, and IV (in vitro) (PubMed:3198760, PubMed:2033050, PubMed:28240246). By cleaving and activating receptor F2RL1/PAR-2, enhances endothelial cell barrier function and thus vascular integrity during neutrophil transendothelial migration (PubMed:23202369). May play a role in neutrophil transendothelial migration, probably when associated with CD177 (PubMed:22266279).[UniProtKB/Swiss-Prot Function]

**Species/Host:** HEK293

**Molecular Weight:** The protein has a predicted molecular mass of 26.2 kDa after removal of the signal peptide. The apparent molecular mass of PRTN3-His is approximately 25–35 kDa due to glycosylation.

**Molecular Characterization:** PRTN3(Ala26–Pro256) 6×His tag

**纯化:** The purity of the protein is greater than 85% as determined by SDS-PAGE and Coomassie blue staining.

**Formulation & Reconstitution:** Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% – 8% trehalose is added as protectants before lyophilization.

**储存和运输:** Store at –20°C to –80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at –80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.



Figure 1. Human PRTN3 Protein, His Tag on SDS-PAGE under reducing condition.