

HUMAN RHOB PROTEIN, HFC TAG

目录: 11901

产品名称: Human RHOB Protein

规格: 10 µg, 50 µg and 100 µg

基因符号: ARH6;ARHB;MST081;MSTP081;RHOH6

Target: RHOB

UNIPROT ID: P62745

描述: Recombinant Human RHOB Protein with C-terminal human Fc tag

背景: Mediates apoptosis in neoplastically transformed cells after DNA damage. Not essential for development but affects cell adhesion and growth factor signaling in transformed cells. Plays a negative role in tumorigenesis as deletion causes tumor formation. Involved in intracellular protein trafficking of a number of proteins. Targets PKN1 to endosomes and is involved in trafficking of the EGF receptor from late endosomes to lysosomes. Also required for stability and nuclear trafficking of AKT1/AKT which promotes endothelial cell survival during vascular development. Serves as a microtubule-dependent signal that is required for the myosin contractile ring formation during cell cycle cytokinesis. Required for genotoxic stress-induced cell death in breast cancer cells.[UniProtKB/Swiss-Prot Function]

Species/Host: HEK293

Molecular Weight: The protein has a predicted molecular mass of 47.4 kDa after removal of the signal peptide. The apparent molecular mass of RHOB-hFc is approximately 35-55 kDa due to glycosylation.

Molecular Characterization: RHOB(Met1-Gly188) hFc(Glu99-Ala330)

纯化: The purity of the protein is greater than 95% as determined by SDS-PAGE and Coomassie blue staining.

Formulation & Reconstitution: Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% – 8% trehalose is added as protectants before lyophilization.

储存和运输: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.



Figure 1. Human RHOB Protein, hFc Tag on SDS-PAGE under reducing condition.