

**HUMAN SLPI PROTEIN, HFC TAG**

目录: 11584

产品名称: Human SLPI Protein

规格: 10 µg, 50 µg and 100 µg

基因符号: ALP;BLPI;HUSI-1;MPI;WAP4;WFDC4

**Target:** SLPI

**UNIPROT ID:** P03973

**描述:** Recombinant human SLPI protein with C-terminal human Fc tag

**背景:** This gene encodes a secreted inhibitor which protects epithelial tissues from serine proteases. It is found in various secretions including seminal plasma, cervical mucus, and bronchial secretions, and has affinity for trypsin, leukocyte elastase, and cathepsin G. Its inhibitory effect contributes to the immune response by protecting epithelial surfaces from attack by endogenous proteolytic enzymes. This antimicrobial protein has antibacterial, antifungal and antiviral activity. [provided by RefSeq, Nov 2014]

**Species/Host:** HEK293

**Molecular Weight:** The protein has a predicted molecular mass of 37.9 kDa after removal of the signal peptide. The apparent molecular mass of SLPI-hFc is approximately 35-55 kDa due to glycosylation.

**Molecular Characterization:** SLPI(Ser26-Ala132) hFc(Glu99-Ala330)

**纯化::** The purity of the protein is greater than 95% as determined by SDS-PAGE and Coomassie blue staining.

**Formulation & Reconstitution:** Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% – 8% trehalose is added as protectants before lyophilization.

**储存和运输:** Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.




Figure 1. Human SLPI Protein, hFc Tag on SDS-PAGE under reducing condition.