

## CD138 MAB (INDATUXIMAB)

目录: 28023

产品名称: CD138 (Indatuximab) Monoclonal Antibody

基因符号: SDC1;Syndecan-1;CD138;SYND1;SDC

描述: CD138 (indatuximab) Monoclonal Antibody

背景: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

经过测试的应用: Flow Cyt

推荐稀释度: Flow Cyt 1:100

种属反应性: Chimeric

亚型: IgG4

纯化: Purified from cell culture supernatant by affinity chromatography

种属反应性: Human CD138

成分: Lyophilized from sterile PBS, pH 7.4. 5 % – 8% trehalose is added as protectants before lyophilization.

储存和运输: Store at  $-20^{\circ}\text{C}$  to  $-80^{\circ}\text{C}$  for 12 months in lyophilized form.

After reconstitution, if not intended for use within a month, aliquot and store at  $-80^{\circ}\text{C}$  (Avoid repeated freezing and thawing).




Figure 1. HEK293 cell line transfected with irrelevant protein (B) and human CD138 (A) were surface stained with CD138 neutralizing antibody  $1\mu\text{g}/\text{ml}$  (indatuximab) followed by Alexa 488-conjugated human IgG secondary antibody.

