

TSHR MAB (K1-70 MAB)

目录: 28083

产品名称: TSHR(K1-70)Monoclonal Antibody

基因符号: CHNG1;hTSHR-I;LGR3

描述: TSHR(K1-70) Monoclonal Antibody

背景: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

经过测试的应用: Flow Cyt

推荐稀释度: Flow Cyt 1:100

种属反应性: Homo sapiens

亚型: IgG1

纯化: Purified from cell culture supernatant by affinity chromatography

种属反应性: Human TSHR

成分: Lyophilized from sterile PBS, pH 7.4. 5 % – 8% trehalose is added as protectants before lyophilization.

储存和运输: Store at -20°C to -80°C for 12 months in lyophilized form.

After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).




Figure 1. Flow cytometry analysis with TSHR(K1-70) mAb 15 μ g/mL on Expi293 cells transfected with Human TSHR (Blue histogram) or Expi293 transfected with irrelevant protein (Red histogram).