

KIT(P577-D579 DEL)

KIT(P577-D579 DEL)

货号: 26288

基因符号: C-Kit; CD117; KIT; PBT; SCFR

描述: KIT(P577-D579 DEL) 小鼠单克隆抗体

背景: Tyrosine-protein kinase Kit, also known as proto-oncogene c-Kit or CD117 or Mast/stem cell growth factor receptor (SCFR), is a receptor tyrosine kinase protein that in humans is encoded by the KIT gene. Activating mutations in this gene are associated with gastrointestinal stromal tumors, testicular seminoma, mast cell disease, melanoma, acute myeloid leukemia, while inactivating mutations are associated with the genetic defect piebaldism.

免疫原: A synthetic peptide from the internal region of KIT which includes the deletion of P577-D579, human origin.

经过测试的应用: ELISA, WB, IHC

推荐稀释度:

ELISA: 1:1000-1:5000

WB: 1:500-1:1000

IHC: 1:50-1:100

浓度: 1 mg/ml

种属反应性: Mouse

形式: Liquid

克隆性: Monoclonal

亚型: IgG

纯化: Purified from ascites


Preservative: No

成分: PBS (without Mg²⁺ and Ca²⁺), pH 7.4, 150 mM NaCl, 50% glycerol

种属反应性: Recognizes P577-D579 DEL mutant, but not wild type KIT of vertebrates.

储存条件: Store at -20°C. Avoid repeated freezing and thawing

Western blot:




1: Kit P577-D579del protein
2: Kit WT protein
WB: anti-Ras (Q61R) mAb

Western blot analysis of recombinant KIT(P577-D579 DEL) and wild type proteins.

Purified His-tagged KIT(P577-D579 DEL) protein (lane 1) and corresponding wild type protein (lane 2) were blotted with KIT(P577-D579 DEL) monoclonal antibody (货号26288).

Immunofluorescence:



Immunofluorescence of cells expressing KIT proteins with KIT(P577-D579 DEL) antibody. HEK293T cells were transfected with pCDNA3-GFP-KIT (WT) plasmid (left column) or pCDNA3-GFP-KIT(P577-D579 DEL) plasmid (right column), then fixed and stained with KIT(P577-D579 DEL) monoclonal antibody (货号26288).