

LGALS1 (DMC277) IGG1 CHIMERIC MAB

目录: 28185

产品名称: LGALS1(DMC277) IgG1 Chimeric Monoclonal Antibody

基因符号: GAL1; GBP

描述: LGALS1 antibody(DMC277) IgG1 Chimeric Monoclonal Antibody

背景: The galectins are a family of beta-galactoside-binding proteins implicated in modulating cell-cell and cell-matrix interactions. This gene product may act as an autocrine negative growth factor that regulates cell proliferation.

经过测试的应用: Flow Cyt

推荐稀释度: Flow Cyt 1:100

种属反应性: Rabbit

亚型: Rabbit:Human Fc chimeric IgG1

纯化: Purified from cell culture supernatant by affinity chromatography

种属反应性: Human LGALS1

成分: Lyophilized from sterile PBS, pH 7.4. 5 % – 8% trehalose is added as protectants before lyophilization.

储存和运输: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).




Figure 1. Flow cytometry analysis with LGALS1 (DMC277) on Expi293 cells transfected with human LGALS1 (Blue histogram) or Expi293 transfected with irrelevant protein (Red histogram).