

## IL7RA MAB (LUSVERTIKIMAB)

目录: 28087

产品名称: IL7RA(Lusvertikimab) Monoclonal Antibody

基因符号: CD127;CDW127;IL-7R- $\alpha$ ;IL7RA;ILRA

描述: IL7RA(Lusvertikimab) Monoclonal Antibody

背景: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

经过测试的应用: Flow Cyt

推荐稀释度: Flow Cyt 1:100

种属反应性: Humanized

亚型: IgG4

纯化: Purified from cell culture supernatant by affinity chromatography

种属反应性: Human IL7RA

成分: Lyophilized from sterile PBS, pH 7.4. 5 % – 8% trehalose is added as protectants before lyophilization.

储存和运输: Store at  $-20^{\circ}\text{C}$  to  $-80^{\circ}\text{C}$  for 12 months in lyophilized form.

After reconstitution, if not intended for use within a month, aliquot and store at  $-80^{\circ}\text{C}$  (Avoid repeated freezing and thawing).




Figure 1. Flow cytometry analysis with 1  $\mu\text{g}/\text{mL}$  IL7RA (Lusvertikimab biosimilar) mAb (28087) on Expi293 cells transfected with Human IL7RA protein (Blue histogram) or Expi293 transfected with irrelevant protein (Red histogram).

