

IL7RA MAB (LUSVERTIKIMAB)

Cat.#: 28087

Product Name: Anti-IL7RA(Lusvertikimab) Monoclonal Antibody

Synonyms: CD127;CDW127;IL-7R-alpha;IL7RA;ILRA

Description: Anti-IL7RA(Lusvertikimab) Monoclonal Antibody

Background: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

Applications: Flow Cyt

Recommended Dilutions: Flow Cyt 1:100

Host Species: Humanized

Isotype: IgG4

Purification: Purified from cell culture supernatant by affinity chromatography

Species Reactivity: Human IL7RA

Constituents: Lyophilized from sterile PBS, pH 7.4. 5 % – 8% trehalose is added as protectants before lyophilization.

Storage & Shipping: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).




Figure 1. Flow cytometry analysis with 1 $\mu\text{g}/\text{mL}$ Anti-IL7RA (lusvertikimab biosimilar) mAb (28087) on Expi293 cells transfected with Human IL7RA protein (Blue histogram) or Expi293 transfected with irrelevant protein (Red histogram).



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
