

ANTI-MCHERRY NANOBODY

mCherry Nanobody

货号: 26401

规格: 100 µL

基因符号:mCherry

描述: mCherry Nanobody

背景: mCherry is a red fluorescent dye widely used in biotechnology as a tracer, including labeling of molecules and localization of cellular components. In biochemical analysis methods such as mass spectrometry and enzyme activity measurement, proteins fused to mCherry and their interacting factors can be quickly and efficiently separated by combining with mCherry-Nanobody.

经过测试的应用: FACS, Elisa, IP, Co-IP, IHC

Dilution: 1:1000

Molecular Weight: 15kD

种属反应性: Human

形式: Liquid

克隆性: Monoclonal

Tag: His, Myc

纯化:: Purified from ascites

Preservative: No

成分: PBS, 40% Glycerol, 0.03% proclin300

Species: Human

储存条件:: Store at 4℃, 6 months; -20℃, 2 years

ELISA mCherry antibody at 1/1000 dilution: 