

## MOUSE 5T4 PROTEIN, HIS TAG

目录: 12211

产品名称: Mouse 5T4 Protein

规格: 10 µg, 50 µg and 100 µg

基因符号: TPBG;WAIF1

**Target:** 5T4

**UNIPROT ID:** Q9Z0L0

**描述:** Recombinant mouse 5T4 protein with C-terminal 6xHis tag

**背景:** Acts upstream of or within several processes, including dendrite arborization; olfactory learning; and positive regulation of intracellular signal transduction. Located in several cellular components, including axon terminus; dendrite; and endoplasmic reticulum. Is expressed in several structures, including brain; epithelium; extraembryonic component; skin; and urinary system. Orthologous to human TPBG (trophoblast glycoprotein). [provided by Alliance of Genome Resources, Apr 2022]

**Species/Host:** HEK293

**Molecular Weight:** The protein has a predicted molecular mass of 36.6 kDa after removal of the signal peptide.

**Molecular Characterization:** Mouse 5T4(Ser32-Ser361) 6×His tag

**纯化:** The purity of the protein is greater than 85% as determined by SDS-PAGE and Coomassie blue staining.

**Formulation & Reconstitution:** Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% – 8% trehalose is added as protectants before lyophilization.

**储存和运输:** Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.




Figure 1. Mouse 5T4 Protein, His Tag on SDS-PAGE under reducing condition.