

**MOUSE ADAM9 PROTEIN, HIS TAG**

目录: 12130

产品名称: Mouse ADAM9 Protein

规格: 10 µg, 50 µg and 100 µg

基因符号: AU020942;MDC9;mKIAA0021;Mltng

**Target:** ADAM9

**UNIPROT ID:** Q61072

**描述:** Recombinant mouse ADAM9 protein with C-terminal 6xHis tag

**背景:** Cleaves and releases a number of molecules with important roles in tumorigenesis and angiogenesis, such as TEK, KDR, EPHB4, CD40, VCAM1 and CDH5 (PubMed:19273593). May mediate cell-cell, cell-matrix interactions and regulate the motility of cells via interactions with integrins (PubMed:10825303).[UniProtKB/Swiss-Prot Function]

**Species/Host:** HEK293

**Molecular Weight:** The protein has a predicted molecular mass of 74.3 kDa after removal of the signal peptide. The apparent molecular mass of mADAM9-His is approximately 55–70 kDa due to glycosylation.

**Molecular Characterization:** Mouse ADAM9(Gly30–Asp697) 6×His tag

**纯化:** The purity of the protein is greater than 85% as determined by SDS-PAGE and Coomassie blue staining.

**Formulation & Reconstitution:** Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% – 8% trehalose is added as protectants before lyophilization.

**储存和运输:** Store at –20°C to –80°C for 12 months in lyophilized form.

After reconstitution, if not intended for use within a month, aliquot and store at –80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.



Figure 1. Mouse ADAM9 Protein, His Tag on SDS-PAGE under reducing condition.