

MOUSE AFP PROTEIN, HIS TAG

目录: 12214

产品名称: Mouse AFP Protein

规格: 10 µg, 50 µg and 100 µg

基因符号: Alpha-fetoprotein

Target: AFP

UNIPROT ID: P02772

描述: Recombinant mouse AFP protein with C-terminal 6xHis tag

背景: Predicted to enable fatty acid binding activity and zinc ion binding activity. Acts upstream of or within SMAD protein signal transduction; ovulation from ovarian follicle; and progesterone metabolic process. Located in cytoplasm and extracellular space. Is expressed in several structures, including amniotic fluid; egg cylinder; gut; liver and biliary system; and metanephros. Human ortholog(s) of this gene implicated in hepatitis C. Orthologous to human AFP (alpha fetoprotein). [provided by Alliance of Genome Resources, Apr 2022]

Species/Host: HEK293

Molecular Weight: The protein has a predicted molecular mass of 66.0 kDa after removal of the signal peptide. The apparent molecular mass of mAFP-His is approximately 55-70 kDa due to glycosylation.

Molecular Characterization: Mouse AFP(Leu21-Val605) 6×His tag

纯化:: The purity of the protein is greater than 85% as determined by SDS-PAGE and Coomassie blue staining.

Formulation & Reconstitution: Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% – 8% trehalose is added as protectants before lyophilization.

储存和运输: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.



Figure 1. Mouse AFP Protein, His Tag on SDS-PAGE under reducing condition.