

## MOUSE BCMA PROTEIN, HFC TAG

目录: 12125

产品名称: Mouse BCMA Protein

规格: 10 µg, 50 µg and 100 µg

基因符号: BCM;BCMA;Tnfrsf13;Tnfrsf13a

**Target:** BCMA

**UNIPROT ID:** O88472

**描述:** Recombinant mouse BCMA protein with C-terminal human Fc tag

**背景:** Receptor for TNFSF13B/BLYS/BAFF and TNFSF13/APRIL. Promotes B-cell survival and plays a role in the regulation of humoral immunity. Activates NF-kappa-B and JNK (By similarity).[UniProtKB/Swiss-Prot Function]

**Species/Host:** HEK293

**Molecular Weight:** The protein has a predicted molecular mass of 31.6 kDa after removal of the signal peptide. The apparent molecular mass of mBCMA-hFc is approximately 35–55 kDa due to glycosylation.

**Molecular Characterization:** Mouse BCMA(Met1–Thr49) hFc(Glu99–Ala330)

**纯化:** The purity of the protein is greater than 95% as determined by SDS-PAGE and Coomassie blue staining.

**Formulation & Reconstitution:** Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% – 8% trehalose is added as protectants before lyophilization.

**储存和运输:** Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.



Figure 1. Mouse BCMA Protein, hFc Tag on SDS-PAGE under reducing condition.




Figure 2. ELISA plate pre-coated by 1 µg/ml (100 µl/well) Mouse BCMA Protein, hFc Tag (12125) can bind anti BCMA monoclonal antibody 28288 in a linear range of 3.2–400 ng/mL.

