

MOUSE PD-1 PROTEIN, HFC TAG

目录: 12221

产品名称: Mouse PD-1 Protein

规格: 10 µg, 50 µg and 100 µg

基因符号: Pdccl1;Pdcl;Ly101

Target: PD-1

UNIPROT ID: Q02242

描述: Recombinant mouse PD-1 protein with C-terminal human Fc tag

背景: Involved in negative regulation of immune response. Acts upstream of or within negative regulation of apoptotic process; negative regulation of tolerance induction; and positive regulation of apoptotic process. Located in external side of plasma membrane. Is expressed in retina. Used to study dilated cardiomyopathy and systemic lupus erythematosus. Human ortholog(s) of this gene implicated in autoimmune disease (multiple); hepatitis B; hepatitis C; hepatocellular carcinoma; and lupus nephritis. Orthologous to human PDCDI (programmed cell death 1). [provided by Alliance of Genome Resources, Apr 2022]

Species/Host: HEK293

Molecular Weight: The protein has a predicted molecular mass of 42.3 kDa after removal of the signal peptide. The apparent molecular mass of mPD-1-hFc is approximately 35–70 kDa due to glycosylation.

Molecular Characterization: Mouse PD-1(Leu25–Gln167) hFc(Glu99–Ala330)

纯化: The purity of the protein is greater than 95% as determined by SDS-PAGE and Coomassie blue staining.

Formulation & Reconstitution: Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% – 8% trehalose is added as protectants before lyophilization.

储存和运输: Store at –20°C to –80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at –80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.



Figure 1. Mouse PD-1 Protein, hFc Tag on SDS-PAGE under reducing condition.