

## MOUSE RNASE4 PROTEIN, HIS TAG

**Cat.#:** 12156

**Product Name:** Mouse RNASE4 Protein

**Size:** 10 µg, 50 µg and 100 µg

**Synonyms:** RNase 4;Rnase4;Ribonuclease 4

**Target:** RNASE4

**UNIPROT ID:** Q9JJH1

**Description:** Recombinant mouse RNASE4 protein with C-terminal 6xHis tag

**Background:** This gene encodes a member of the pancreatic ribonuclease A superfamily. The encoded enzyme is secreted and has unique uridine specificity. This gene resides in a cluster of highly related genes. It shares dual promoters and 5' exons with the angiogenin, ribonuclease, RNase A family, 5 gene. Each gene splices to a unique downstream exon that contains its complete coding region. Two alternatively spliced variants, with different 5' exons but the same coding exon, have been identified. [provided by RefSeq, Jun 2009]

**Species/Host:** HEK293

**Molecular Weight:** The protein has a predicted molecular mass of 14.7 kDa after removal of the signal peptide. The apparent molecular mass of mRNASE4-His is approximately 10–15 kDa due to glycosylation.

**Molecular Characterization:** Mouse RNASE4(Gln30-Cys472) 6xHis tag

**Purity:** The purity of the protein is greater than 85% as determined by SDS-PAGE and Coomassie blue staining.

**Formulation & Reconstitution:** Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% – 8% trehalose is added as protectants before lyophilization.

**Storage & Shipping:** Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.




Figure 1. Mouse RNASE4 Protein, His Tag on SDS-PAGE under reducing condition.