

MOUSE VWF (764-2813) PROTEIN, HIS TAG

目录: 12149

产品名称: Mouse VWF (764-2813) Protein

规格: 10 µg, 50 µg and 100 µg

基因符号: F8VWF;VWD

Target: VWF

UNIPROT ID: Q8CIZ8

描述: Recombinant mouse VWF(764-2813) protein with C-terminal 6xHis tag

背景: This gene encodes a glycoprotein involved in hemostasis. The encoded preproprotein is proteolytically processed following assembly into large multimeric complexes. These complexes function in the adhesion of platelets to sites of vascular injury and the transport of various proteins in the blood. Mutations in this gene result in von Willebrand disease, an inherited bleeding disorder. An unprocessed pseudogene has been found on chromosome 22. [provided by RefSeq, Oct 2015]

Species/Host: HEK293

Molecular Weight: The protein has a predicted molecular mass of 225.9 kDa after removal of the signal peptide.

Molecular Characterization: Mouse VWF(Ser764-Lys2813) 6×His tag

纯化: The purity of the protein is greater than 85% as determined by SDS-PAGE and Coomassie blue staining.

Formulation & Reconstitution: Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% – 8% trehalose is added as protectants before lyophilization.

储存和运输: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.




Figure 1. Mouse VWF (764-2813) Protein, His Tag on SDS-PAGE under reducing condition.