

## CD37 MAB (NARATUXIMAB)

目录: 28055

产品名称: CD37 (Naratuximab) Monoclonal Antibody

基因符号: CD37;TSPAN26;Tspan-26

描述: CD37 (naratuximab ) Monoclonal Antibody

背景: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

经过测试的应用: Flow Cyt

推荐稀释度: Flow Cyt 1:100

种属反应性: Chimeric

亚型: IgG1

纯化: Purified from cell culture supernatant by affinity chromatography

种属反应性: Human CD37

成分: Lyophilized from sterile PBS, pH 7.4. 5 % – 8% trehalose is added as protectants before lyophilization.

储存和运输: Store at -20°C to -80°C for 12 months in lyophilized form.

After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).




Figure 1. HEK293 cell line transfected with irrelevant protein (B) and human CD37 (A) were surface stained with CD37 mAb 1 µg/ml (naratuximab) followed by Alexa 488-conjugated human IgG secondary antibody.

