

## MASP2 MAB (NARSOPLIMAB)

目录: 28135

产品名称: MASP2(Narsoplimab) Monoclonal Antibody

基因符号: MASP-2

描述: MASP2(narsoplimab) Monoclonal Antibody

背景: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

经过测试的应用: Flow Cyt

推荐稀释度: Flow Cyt 1:100

种属反应性: Homo sapiens

亚型: IgG4

纯化: Purified from cell culture supernatant by affinity chromatography

种属反应性: Human MASP2

成分: Lyophilized from sterile PBS, pH 7.4. 5 % – 8% trehalose is added as protectants before lyophilization.

储存和运输: Store at -20°C to -80°C for 12 months in lyophilized form.

After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).




Figure 1. MASP2 protein is highly expressed on the surface of Jurkat cell membrane. Flow cytometry analysis with 15 µg/mL MASP2(narsoplimab biosimilar) mAb (28135) (Blue histogram) or isotype control mAb (Red histogram) on Jurkat cells.

