

HUMAN OR52D1 FULL LENGTH PROTEIN

目录: 11118

产品名称: Human OR52D1 Full Length Protein

规格: 10 µg; 50 µg and 100 µg

基因符号: OR11-43

Target: OR52D1

UNIPROT ID: Q9H346

描述: Human OR52D1 Full Length Protein-Synthetic Nanodisc

背景: Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms.

Species/Host: HEK293

Molecular Weight: The human full length OR52D1 protein has a MW of 35.1 kDa

Formulation & Reconstitution: Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% - 8% trehalose is added as protectants before lyophilization.

储存和运输: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.

Protein Families: Transmembrane

Protein Pathways: Olfactory transduction

ELISA assay to evaluate OR52D1-Nanodisc 0.2µg Human OR52D1-Nanodisc per well




Figure 1. Elisa plates were pre-coated with Flag Tag OR52D1-Nanodisc (0.2µg/per well). Serial diluted Flag monoclonal antibody solutions were added, washed, and incubated with secondary antibody before Elisa reading. From above data, the EC50 for Flag monoclonal antibody binding with OR52D1-Nanodisc is 3.737ng/ml.




Figure 2. WB analysis of Human OR52D1-Nanodisc with Flag monoclonal antibody at 1/5000 dilution, followed by Goat Rabbit IgG HRP at 1/5000 dilution