

TFRC MAB (PABINAFUSP ALFA MAB)

Cat.#: 28132

Product Name: Anti-TFRC(Pabinafusp Alfa) Monoclonal Antibody

Synonyms: TR;TfR;TfR1;Trfr;T9;p90;CD71

Description: Anti-TFRC(pabinafusp alfa) Monoclonal Antibody

Background: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

Applications: Flow Cyt

Recommended Dilutions: Flow Cyt 1:100

Host Species: Humanized

Isotype: IgG1

Purification: Purified from cell culture supernatant by affinity chromatography

Species Reactivity: Human TFRC

Constituents: Lyophilized from sterile PBS, pH 7.4. 5 % – 8% trehalose is added as protectants before lyophilization.

Storage & Shipping: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).




Figure 1. Flow cytometry analysis with 1 $\mu\text{g}/\text{mL}$ Anti-TFRC mAb (28132) on Expi293 cells transfected with Human TFRC protein (Blue histogram) or Expi293 transfected with irrelevant protein (Red histogram).