

## **TFRC MAB (PABINAFUSP ALFA MAB)**

目录: 28132

产品名称: TFRC(Pabinafusp Alfa) Monoclonal Antibody

基因符号: TR;TfR;TfR1;Trfr;T9;p90;CD71

描述: TFRC(pabinafusp alfa ) Monoclonal Antibody

背景: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

经过测试的应用: Flow Cyt

推荐稀释度: Flow Cyt 1:100

种属反应性: Humanized

亚型: IgG1

纯化: Purified from cell culture supernatant by affinity chromatography

种属反应性: Human TFRC

成分: Lyophilized from sterile PBS, pH 7.4. 5 % – 8% trehalose is added as protectants before lyophilization.

储存和运输: Store at -20°C to -80°C for 12 months in lyophilized form.

After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).




Figure 1. Flow cytometry analysis with 1  $\mu$ g/mL TFRC mAb (28132) on Expi293 cells transfected with Human TFRC protein (Blue histogram) or Expi293 transfected with irrelevant protein (Red histogram).