

CCN2 MAB (PAMREVLUMAB)

目录: 28105

产品名称: CCN2(Pamrevlumab) Monoclonal Antibody

基因符号: CTGF;HCS24;IGFBP8;NOV2

描述: CCN2(pamrevlumab) Monoclonal Antibody

背景: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

经过测试的应用: Flow Cyt

推荐稀释度: Flow Cyt 1:100

种属反应性: Homo sapiens

亚型: IgG1

纯化: Purified from cell culture supernatant by affinity chromatography

种属反应性: Human CCN2

成分: Lyophilized from sterile PBS, pH 7.4. 5 % – 8% trehalose is added as protectants before lyophilization.

储存和运输: Store at -20°C to -80°C for 12 months in lyophilized form.

After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).




Figure 1. Flow cytometry analysis under cell membrane permeable condition with 1 μ g/mL CCN2 (pamrevlumab biosimilar) mAb (28105) on Expi293 cells transfected with Human CCN2 protein (Blue histogram) or Expi293 transfected with irrelevant protein (Red histogram).