

PE-CONJUGATED ANTI-GPRC5D (DM91) RABBIT MAB

目录: 28425

产品名称: PE-Conjugated GPRC5D(DM91) Rabbit Monoclonal Antibody

基因符号: GPRC5D

描述: PE-conjugated GPRC5D antibody(DM91) Rabbit Monoclonal Antibody

背景: The protein encoded by this gene is a member of the G protein-coupled receptor family; however; the specific function of this gene has not yet been determined.

经过测试的应用: Flow Cyt

推荐稀释度: Flow Cyt 1:100

种属反应性: Rabbit

亚型: Rabbit IgG

纯化: Purified from cell culture supernatant by affinity chromatography

种属反应性: Human GPRC5D

成分: Lyophilized from sterile PBS, pH 7.4. 5 % – 8% trehalose is added as protectants before lyophilization.

储存和运输: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).




Figure 1. Flow cytometry analysis with 100 μ l of 1:100 dilution PE-conjugated GPRC5D(DM91) mAb (28425) on Expi293 cells transfected with Human GPRC5D protein (Blue histogram) or Expi293 transfected with irrelevant protein (Red histogram).

