

## PE-CONJUGATED HUMAN GPRC5D PROTEIN

目录: 12235

产品名称: PE-Conjugated Human GPRC5D Protein

规格: 10 µg, 50 µg and 100 µg

基因符号: GPRC5D

**Target:** GPRC5D

**UNIPROT ID:** Q9NZD1

**描述:** PE-conjugated Recombinant human GPRC5D protein with C-terminal mouse Fc tag

**背景:** The protein encoded by this gene is a member of the G protein-coupled receptor family, however, the specific function of this gene has not yet been determined.

**Species/Host:** HEK293

**Molecular Weight:**

**Molecular Characterization:** GPRC5D(Met2-Glu21) mFc(Pro99-Lys330)

**纯化::** The purity of the protein is greater than 95% as determined by SDS-PAGE and Coomassie blue staining.

**Formulation & Reconstitution:** Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% – 8% trehalose is added as protectants before lyophilization.

**储存和运输:** Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.




Figure 1. Flow cytometry analysis with 1 µg/ml PE-conjugated GPRC5D-mFc protein (12235) on Anti GPRC5D CAR-Jurkat cells (Blue histogram) or negative Jurkat cells (Red histogram).