

PIK3CA(G1049S)

PIK3CA(G1049S)

货号: 26219

规格: 100 μ L

基因符号: PIK3CA, CLOVE, CWS5, MCAP, MCM, MCMTTC, PI3K, p110- α

描述: PIK3CA(G1049S) 小鼠单克隆抗体

背景: Phosphatidyl 3-kinases (PI3K) are a family of lipid kinases involved in many cellular processes, including cell growth, proliferation, differentiation, motility, and survival. It is composed of an 85 kDa regulatory subunit and a 110 kDa catalytic subunit. Mutants of PIK3CA has been found in the pathogenesis of several cancers, including colon cancer, gliomas, gastric cancer, breast cancer, endometrial cancer, and lung cancer.

免疫原: A synthetic peptide from the internal region of PIK3CA which includes the mutation of G1049S, human origin.

经过测试的应用: ELISA, WB, IHC

推荐稀释度:

ELISA: 1:1000-1:5000

WB: 1:500-1:2000

浓度: 1 mg/ml

种属反应性: Mouse

形式: Liquid

克隆性: Monoclonal

亚型: IgG

纯化:: Purified from ascites

Preservative: No

成分: PBS (without Mg^{2+} and Ca^{2+}), pH 7.4, 150 mM NaCl, 50% glycerol


种属反应性: Recognizes G1049S mutant, but not wild type PIK3CA of vertebrates.

储存条件:: Store at $-20^{\circ}C$. Avoid repeated freezing and thawing

Western blot:

PIK3CA protein

WT G1049S



WB: Anti- PIK3CA(G1049S) mAb

Western blot analysis of recombinant PIK3CA(G1049S) and wild type proteins. Purified His-tagged PIK3CA(G1049S) protein (lane 2) and corresponding wild type protein (lane 1) were blotted with PIK3CA(G1049S) monoclonal antibody (货号26219).