

RAB17 PAB

Rab17 Polyclonal Antibody

货号: 21049

基因符号: RAB17

描述: Rab17 Rabbit Polyclonal Antibody

背景: The Rab-subfamily of small GTPases plays an important role in the regulation of membrane trafficking. Rab17 is an epithelial cell-specific GTPase, which is involved in transcellular transport.

免疫原: A synthetic peptide from the C-terminal end of human Rab17 protein

经过测试的应用: ELISA, WB, IHC

推荐稀释度:

ELISA: 1:100-1:5000

WB: 1:100-1:2000

IHC: 1:100-1:1000

浓度: 0.2 mg/ml

种属反应性: Rabbit

形式: Liquid

克隆性: Polyclonal

纯化: Purified from serum


Preservative: No

成分: PBS (without Mg^{2+} and Ca^{2+}), pH 7.4, 150 mM NaCl, 50% glycerol

种属反应性: Rab17 antibody recognizes Rab17 from vertebrates.

储存条件: Store at $-20^{\circ}C$. Avoid repeated freezing and thawing

Western blot:



WB: anti-Rab17 pAb

HepG whole-cell lysate was used for Western blot with anti Rab17 polyclonal antibody (cat # 21049).