

## RAB35 MAB

### RAB35 Monoclonal Antibody

货号: 26033

基因符号: RAB35

描述: Rab35 小鼠单克隆抗体

背景: The Rab family of GTPases regulates eukaryotic vesicular membrane traffic. Rab35 is with both plasma membrane and endosomal localization, and has been implicated in diverse processes that include T-cell receptor recycling, oocyte yolk protein recycling and cytokinesis. Rab35 regulates neurite outgrowth in neuronal-like cells, and can induce protrusions.

免疫原: Recombinant full length protein of Rab35

经过测试的应用: ELISA, WB, IHC, IP

推荐稀释度:

ELISA: 1:100-1:5000

WB: 1:100-1:2000

IHC: 1:100-1:1000

浓度:

种属反应性: Mouse

形式: Liquid

克隆性: Monoclonal

亚型: IgG

纯化:: Purified from ascites


**Preservative:** No

成分: PBS (without  $Mg^{2+}$  and  $Ca^{2+}$ ), pH 7.4, 150 mM NaCl, 50% glycerol

种属反应性: Rab35 antibody recognizes Rab35 of vertebrates.

储存条件:: Store at  $-20^{\circ}C$ . Avoid repeated freezing and thawing

**Western blot:**



**WB:anti-Rab35 mAb**

**Western blot analysis of recombinant GST-Rab35 proteins.**