

RAB4A PAB

Rab4A Polyclonal Antibody

货号: 21158

基因符号: RAB4A

描述: Rab4A Rabbit Polyclonal Antibody

背景: Rab-4A belongs to the small GTPase Rab family. It interacts with RAB11FIP1, RABEP1, ZFYVE20, RUFY1, SGSM1, SGSM2 and SGSM3. It participates in protein transport and in vesicular traffic.

免疫原: Purified Rab4A protein.

经过测试的应用: ELISA, WB, IHC, IP

推荐稀释度:

ELISA: 1:100-1:5000

WB: 1:100-1:2000

IHC: 1:10-1:100

浓度: 0.15 mg/ml

种属反应性: Rabbit

形式: Liquid

克隆性: Polyclonal

亚型: IgG

纯化:: Purified from serum

Preservative: No

成分: PBS (without Mg^{2+} and Ca^{2+}), pH 7.4, 150 mM NaCl, 59% glycerol Species

种属反应性: Rab4A antibody recognizes Rab4A proteins from vertebrates.

储存条件:: Store at $-20^{\circ}C$. Avoid repeated freezing and thawing

Product Description

Immunoprecipitation/Western blot:



Purified GST-tagged Rab4A proteins (货号10130) were used for Western blot with anti Rab4A polyclonal antibody (货号21158).