

RAB4A PAB

Rab4A Polyclonal Antibody

Cat. #: 21158

Gene Symbol: RAB4A

Description: Anti-Rab4A Rabbit Polyclonal Antibody

Background: Rab-4A belongs to the small GTPase Rab family. It interacts with RAB11FIP1, RABEP1, ZFYVE20, RUFY1, SGSM1, SGSM2 and SGSM3. It participates in protein transport and in vesicular traffic.

Immunogen: Purified Rab4A protein.

Applications: ELISA, WB, IHC, IP

Recommended Dilutions:

ELISA: 1:100–1:5000

WB: 1:100–1:2000

IHC: 1:10–1:100

Concentration: 0.15 mg/ml

Host Species: Rabbit

Format: Liquid

Clonality: Polyclonal

Isotype: IgG

Purity: Purified from serum

Preservative: No

Constituents: PBS (without Mg^{2+} and Ca^{2+}), pH 7.4, 150 mM NaCl, 59% glycerol Species

Species Reactivity: Rab4A antibody recognizes Rab4A proteins from vertebrates.

Storage Conditions: Store at $-20^{\circ}C$. Avoid repeated freezing and thawing

Immunoprecipitation/Western blot:



Purified GST-tagged Rab4A proteins (Cat. #10130) were used for Western blot with anti Rab4A polyclonal antibody (Cat. #21158).