

RAB5 MAB

Rab5

货号 26016

基因符号: RAB5

描述: Rab5 小鼠单克隆抗体

背景: Rab5 is a member of Rab-family GTPases. Rab5 is a rate-limiting component of the machinery regulating the kinetics of membrane traffic in the early endocytic pathway.

免疫原: Recombinant full length protein of Rab5A

经过测试的应用: ELISA, WB, IHC

推荐稀释度:

ELISA: 1:100-1:5000

WB: 1:100-1:2000

IHC: 1:50-1:100

浓度: 1 mg/ml

种属反应性: mouse

形式: Liquid

克隆性: Monoclonal

亚型: IgG

纯化:: Purified from ascites


Preservative: no

成分: PBS (without Mg^{2+} and Ca^{2+}), pH 7.4, 150 mM NaCl, 50% glycerol

种属反应性: Rab5 antibody recognizes Rab5 proteins of vertebrates.

储存条件:: Store at $-20^{\circ}C$. Avoid repeated freezing and thawing

Western blot:



WB: anti-Rab5 MAb

Western blot analysis of Rab5 proteins. EC-1 whole-cell lysate was used for Western blot.