

RAB5B PAB

Rab5B Polyclonal Antibody

货号: 21037

基因符号: RAB5B

描述: Rab5B Rabbit Polyclonal Antibody

背景: : Rab5, consists of three isoforms, Rab5A, Rab5B, and Rab5C, is a small GTPase that is localized to early endosomes. It regulates the fusion between endocytic vesicles and early endosomes, as well as the homotypic fusion between early endosomes.

免疫原: A synthetic peptide from Rab5B

经过测试的应用: ELISA, WB, IHC

推荐稀释度 : ELISA 1:100-1:5000

WB 1:100-1:2000

IHC 1:100-1:1000

浓度: 0.2 mg/ml

种属反应性: Rabbit

形式: Liquid

克隆性 : Polyclonal

亚型: IgG

纯化:: Purified from serum

Preservative: No

成分: PBS (without Mg^{2+} and Ca^{2+}), pH 7.4, 150 mM NaCl, 50% glycerol

种属反应性: Rab5B antibody recognizes Rab5B from vertebrates.

储存条件:: Store at $-20^{\circ}C$. Avoid repeated freezing and thawing

Western blot:

HepG2



WB: anti-Rab5B pAb

Western blot analysis of Rab5B proteins. HepG2 cell lysate was used for Western blot with anti Rab5B polyclonal antibody (cat # 21037).