

GITR:TNFRSF18 MAB (RAGIFILIMAB)

目录: 28072

产品名称: GITR:TNFRSF18(Ragifilimab) Monoclonal Antibody

基因符号: AITR;GITR;TNFRSF18;CD357

描述: GITR:TNFRSF18(ragifilimab) Monoclonal Antibody

背景: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

经过测试的应用: Flow Cyt

推荐稀释度: Flow Cyt 1:100

种属反应性: Humanized

亚型: IgG4

纯化: Purified from cell culture supernatant by affinity chromatography

种属反应性: Human GITR

成分: Lyophilized from sterile PBS, pH 7.4. 5 % – 8% trehalose is added as protectants before lyophilization.

储存和运输: Store at -20°C to -80°C for 12 months in lyophilized form.

After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).




Figure 1. Flow cytometry analysis with GITR/TNFRSF18(ragifilimab) mAb 1 μ g/mL on Expi293 cells transfected with Human GITR (Blue histogram) or Expi293 transfected with irrelevant protein (Red histogram).