

## RALB MAB

### RalB Monoclonal Antibody

货号: 26018

基因符号: RALB

描述: RalB 小鼠单克隆抗体

背景: The Ras-like small G proteins, RalA/B, are important components of Ras signaling pathways, implicated in the initiation and maintenance of tumorigenic transformation, as well as vesicle transport, apoptosis, transcription, cell migration, and cell proliferation.

免疫原: Recombinant full length protein of RalB

经过测试的应用: ELISA, WB, IP, IHC

推荐稀释度: ELISA: 1:1000-1:5000

WB: 1:100-1:2000

IP: 1 µg for 1~2 mg total cellular proteins

IHC: 1:100

浓度: 1 mg/ml

种属反应性: Mouse

形式: Liquid

克隆性: Monoclonal

亚型: IgG

纯化: Purified from ascites

**Preservative:** No

成分: PBS (without  $Mg^{2+}$  and  $Ca^{2+}$ ), pH 7.4, 150 mM NaCl, 50% glycerol

种属反应性: RalB antibody recognizes RalB of vertebrates.

储存条件: Store at -20°C. Avoid repeated freezing and thawing

### Western blot:

**Western blot:**

**kD**

**43 –**

**33 –**

**26 –**



**RalB**

**WB: anti-RalB mAb**

Western blot analysis of RalB proteins. EC-1 cell lysate was used for Western blot with RalB monoclonal antibody (cat. # 26018).