

## RALB PAB

### RalB Polyclonal Antibody

**Cat. #:** 21033

**Gene Symbol:** RALB

**Description:** Anti-RalB Rabbit Polyclonal Antibody

**Background:** The Ras-like small G proteins, RalA/B, are important components of Ras signaling pathways, implicated in the initiation and maintenance of tumorigenic transformation, as well as vesicle transport, apoptosis, transcription, cell migration, and cell proliferation.

**Immunogen:** A synthetic peptide from the internal region of RalB protein.

**Applications:** ELISA, WB, IHC

**Recommended Dilutions:** ELISA 1:100–1:5000

WB 1:100–1:2000

IHC 1:100–1:1000

**Concentration:** 0.2 mg/ml

**Host Species:** Rabbit

**Format:** Liquid

**Clonality:** Polyclonal

**Isotype:** IgG

**Purity:** Purified from serum


**Preservative:** No

**Constituents:** PBS (without  $Mg^{2+}$  and  $Ca^{2+}$ ), pH 7.4, 150 mM NaCl, 50% glycerol

**Species Reactivity:** RalB antibody recognizes RalB from vertebrates.

**Storage Conditions:** Store at  $-20^{\circ}C$ . Avoid repeated freezing and thawing

## Western blot:



WB: anti-RalB pAb

## Western blot analysis of RalB proteins.

EC-1 cell lysate was used for Western blot with anti RalB polyclonal antibody (cat # 21033).