

RAP1B PAB**Rap1B Polyclonal Antibody**

货号: 21156

基因符号: RAP1B

描述: Rap1B Rabbit Polyclonal Antibody

背景: Small GTPases are a super-family of cellular signaling regulators. Rap family members (Rap1A, Rap1B, Rap2A, Rap2B, and Rap2C) are Ras-like GTPases that regulate cell proliferation, differentiation, apoptosis and adhesion mechanisms. Rap1 (Ras-proximate-1) regulates apoptosis through interaction with two MAP3K effectors of the MEK/ERK MAP-kinase pathway, Raf-1 and B-Raf. Rap1 regulates cell adhesion mechanisms through effector molecules such as the integrins and cadherins.

免疫原: Purified His-tagged Rap1B protein.

经过测试的应用: ELISA, WB, IHC, IP

推荐稀释度:

ELISA: 1:100-1:5000

WB: 1:100-1:1000

IHC: 1:100-1:1000

浓度: 0.2 mg/ml

种属反应性: Rabbit

形式: Liquid

克隆性: Polyclonal

亚型: IgG

纯化: Purified from serum


Preservative: No

成分: PBS (without Mg^{2+} and Ca^{2+}), pH 7.4, 150 mM NaCl, 50% glycerol

种属反应性: Rap1B antibody recognizes Rap1B proteins from vertebrates.

储存条件: Store at $-20^{\circ}C$. Avoid repeated freezing and thawing

Immunoprecipitation/Western blot:



Purified His-tagged Rap1B protein (货号10120) was used for western blot with the polyclonal anti Rap1B antibody (货号 21156).